

BRUKER AXS SMART APEXII SINGLE CRYSTAL DIFFRACTOMETER

STANDARD OPERATION PROCEDURE



Dr. Chunhua (Tony) Hu
New York University, Department of Chemistry
May 10, 2010

1. SAFETY

Your safety is the first concern. Each individual who wishes to be authorized to use X-ray equipments at the Facility must first receive training concerning the X-ray radiation hazards associated with the use of the equipment, the function and importance of the equipment's safety devices, the proper operation procedures, and procedures to follow in the even of a suspected radiation exposure.

1.1 Links to On-line Safety Training:

http://www.safety.ncsu.edu/xray_training/

<http://www.nyu.edu/fas/dept/chemistry/wardgroup/DoubleEdgedSword.mov>

1.2 Radiation Dosimeter

Wear your radiation badge and finger dosimeter all the time when you work in the X-ray lab.



To obtain the whole body badge and finger dosimeter, please contact Charles Strom (charlie.strom@nyu.edu, 212-998-88480, 10011 Silver).

1.3 Geiger counter

A Geiger counter is available in the lab for detecting any X-ray leakage.



1.4 Standard Operation Protocol

Always follow this standard operation protocol to operate the machine. If you have any questions, do not hesitate to contact Dr. Chunhua (Tony) Hu (212-9988769, chunhua.hu@nyu.edu).

2. DATA COLLECTION

2.1. Log in the computer with the username of *APEXII User* and the password of *apexii*.

2.2. To start BIS, choose **Start > Programs > Bruker AXS > Administration > BIS** or double-click BIS icon on the Desktop. After initialization, BIS asks for the detector distance as shown in the Fig. 1. Click OK.




Fig. 1 BIS login

2.3. To start APEX2, choose **Start > Programs > Bruker AXS > APEX2** or double-click APEX2 icon on the Desktop. Choose **Instrument > Connect**, use the Host Name of *apexii* to connect the diffractometer.





Fig. 2 Connect the host computer

2.4. To create a new sample, choose **Sample > New...** from the Menu Bar or click the New icon . Give a sample name under the pop-up window (Fig. 3). The sample name should be year (10 stands for 2010), initials of the group advisor (three letters, for example, mdw is Mike D. Ward), and the sequence number followed by an initial of the operator (for example, h stands for Chunhua Hu).




Fig. 3 New Sample

2.5. If VIDEO is not started, start it by double-clicking the VIDEO icon on the desktop.

2.6. In APEX2 program choose **Center Crystal**.




Fig. 4 Center Crystal

2.7. Click **HOME** to home the goniometer. Adjust the screw in the goniometer head and move the crystal to the center of the video crosshair. Click **Spin Phi 180** to check the position. Adjust it if necessary. Click **Spin Phi 90** to rotate the Phi in 90 degrees and move the crystal to the center. Click **Spin Phi 180** to check the position again. Adjust it if necessary. Click **Spin Phi 90** and **Spin Phi 180** a few times more to make sure that the crystal is centered.


2.8. Click  icon in the VIDEO program. Measure the crystal size by mouse clicking. Write down the crystal information in the **Describe** module (Fig. 5).

Fig. 5 Sample Description

2.9. Select **Simple Scans** module. Click **Still**, and **Drive + Scan**. After 10 seconds, an image collected is displayed on the screen. If you are not satisfied with the image (too weak intensity, obviously twinned crystal, or split crystal like

in Fig. 6), find another crystal under microscope and try again. If the crystal diffracts well and gives nice spot distribution, go ahead to collect matrix runs and determine the unit cell.




Fig. 6 Simple Scans

2.10. Select **Evaluate** module. Click **Determine Unit Cell**. Click **Run** to start the automatic unit cell determination. The unit cell parameters are listed in the table afterwards (Fig. 7).




Fig. 7 Determine Unit Cell

2.11. Select **Collect** module. Click **Experiment**. Click **Load Table...** and load the file MIT_5cm_5runs.exp under c:\frames\. Change **Default time** (10 seconds) to *shorter* for the strongly diffracting crystal and *longer* for the weak ones. Click **Validate**. If passed, click **Execute** to collect data.




Fig. 8 Experiment

3. DATA PROCESSING

3.1. After data collection, select **Integrate/Integrate Images**.




Fig. 9 Integrate Images

Click **Refinement Options** button. In the pop-up window, enable the options shown in the Fig. 10. Click **OK**.




Fig. 10 Refinement Options

Click **Integrate Options** button. In the pop-up window, click **More Options** to extend the window. Enable the options shown in the Fig. 11 and use 0.500 for the **Generate Mask** option. Click **OK**.

Fig. 11 Integrate Options

Click **Start Integration**. During the integration, the XYZ box size is re-calculated and updated. Write down these numbers, i.e. 0.73 for X, 0.62 for Y, and 1.38 for Z in this example (Fig. 12). When the integration finishes, select **Sample > import > P4P/SPIN file**, load the *.p4p file and import all the information from this file (Fig. 13). Change the **Initial XYZ Box Size** with the updated ones under **Reflection Options**, re-run the integration. Doing so usually gives better results.




Fig. 12 Integration in Progress




Fig. 13 Import P4P file

3.2. Select **Scale** module and click **scale** icon. In most cases the default settings are good enough (Fig. 14). If the space group of the structure is non-centrosymmetric, you need to select the corresponding non-centrosymmetric point group. Click **Next** to open the **Parameter Refinement** window (Fig. 15). Select the proper **Absorption Type** for your crystal. Click **Refine**. A graph (Fig. 16) appeared shows the intensity statistics. Click **Next** to continue. In the **Error Model** window (Fig. 17) click **Determine Error Model**. If you are happy with the results, click **Finish**. In the next window (Fig. 18) click on the tabs at the bottom of the screen to view diagnostic data. If the results are good, click **Exit AXScale** to exit the **Scale**.

Scale

Setup Numerical Absorption Correction Parameter Refinement Error Model Diagnostics

Use Merged Batches or Individual Batches

☒ Merged Batches ☐ Individual Batches

E:\frames\guest\10ktk1h\work\10ktk1h_0m.raw

Base: 10ktk1h

Output File Type: Unmerged .hkl file

Output File Name: 10ktk1h_0m

Diagnostic Plots File Name: 10ktk1h.eps

Title of Diagnostic Plots: 10ktk1h

Log File: 10ktk1h.abs

☒ Use only centrosymmetric point groups

Point Group: 2/m

☒ Additional Spherical Absorption Correction

Mu*r of Equivalent Sphere: 0.2

Allow for crystal decomposition by B-value refinement: None

Extra Linear Correction to be Applied to Each Reflection: None

Spatial display of $|I - \langle I \rangle|/su$ greater than: 3.0

Absorption Correction Type

☒ Multiscan Absorption Correction

☐ Numerical Absorption Correction (from Face Indices)

Next

Finish

Fig. 14 Scale Setup




Fig. 15 Scale Parameter Refinement




Fig. 16 Scale Graph




Fig. 17 Scale Error Model

Scale

Setup Numerical Absorption Correction Parameter Refinement Error Model Diagnostics

Direction cosine

Mean error 0.000

Maximum error 0.002

Data Statistics

Maximum 2-Theta (degrees) 56.71

Maximum Resolution (Angstroms) 0.75

Approximate Wavelength (Angstroms) 0.71076

Reflections 32138

Unique 4185

Data per frame 12.28

Average Redundancy 7.68

Approximate Unit Cell (from direction cosines)

| a | b | c | α | β | γ |
|-------|--------|--------|----------|---------|----------|
| 9.076 | 10.161 | 17.587 | 90.006 | 101.040 | 89.974 |

Unit Cell Parameters From Database

| a | b | c | α | β | γ |
|-------|--------|--------|----------|---------|----------|
| 9.069 | 10.155 | 17.581 | 90.000 | 101.023 | 90.000 |

Initial Reflections

Total 32138

Unique 4185

Error Model g-value

Suggested 0.0478

Use 0.0478

Transmission Data

Corrected Reflections: 32128

Minimum Transmission 0.5932

Maximum Transmission 0.7457

Ratio of min/max apparent transmission 0.7955

wR2(int)

Initial wR2(int) 0.0771

Overall wR2(int) 0.0549

(selected reflections only, after parameter refinement)

Reflections after Outlier Rejection

Total 32128

% Rejected 0.0

Unique 4185

% Rejected 0.0

Statistics Reflection Graphs Refinement Graph Error Model Graphs Scale Variations Intensity Statistics Chi-Squared Spatial Distribution Spatial Distribution 2 Spatial Distribution 3

Start Over Exit AXXScale

Fig. 18 Scale Diagnostics

3.3 Use the obtained *.p4p and *.hkl files as input files for space group determination (XPREP), structure solution (SHELTL xs), and structure refinement (SHELTL xl).

Please refer to the **APEX2 User Manual** for more details.

4. EXPERIMENTAL DESCRIPTION

A needle crystal with the size of 0.05 × 0.05 × 0.40 mm³ was selected for geometry and intensity data collection with a Bruker SMART APEXII CCD area detector on a D8 goniometer at 100 K. The temperature during the data collection was controlled with an Oxford Cryosystems Series 700 plus instrument. Preliminary lattice parameters and orientation matrices were obtained from three sets of frames. Data were collected using graphite-monochromated and 0.5 mm-MonoCap-collimated Mo-K_α radiation ($\lambda = 0.71073$ Å) with the ω scan method [1]. Data were processed with the INTEGRATE program of the APEX2 software [1] for reduction and cell refinement. Multi-scan absorption corrections were applied by using the SCALE program for area detector. The crystal faces were determined by APEX2/Crystal Faces [1]. The structure was solved by the direct method and refined on F² (SHELXTL) [2]. Non-hydrogen atoms were refined with anisotropic displacement parameters, and hydrogen atoms on carbons were placed in idealized positions (C-H = 0.93 or 0.96 Å) and included as riding with $U_{\text{iso}}(\text{H}) = 1.2$ or $1.5 U_{\text{eq}}(\text{non-H})$. Crystal data, data collection parameters and refinement results are listed in Table 1. Crystallographic data (excluding structure factors) for the structure reported in this paper have been deposited with the Cambridge Crystallographic Data Centre as supplementary publications no. ## (**compound 1**). Copies of available material can be obtained, free of charge, on application to CCDC, 12 Union Road, Cambridge CB2 1EZ, UK, (fax: +44-(0)1223-336033 or e-mail: deposit@ccdc.cam.ac.uk).

[1] APEX2 (version 2009.11-0). *Program for Bruker CCD X-ray Diffractometer Control*, Bruker AXS Inc., Madison, WI, 2009.

[2] G. M. Sheldrick, SHELXTL, version 6.14. Program for solution and refinement of crystal structures, Universität Göttingen, Germany, 2000.

The underlined contents need be revised according to the experimental records and the journal requirements. The crystal size and color and the data collection temperature can be found in the CIF:

| | |
|------------------------------------|--------|
| <u>_exptl_crystal_description</u> | block |
| <u>_exptl_crystal_colour</u> | yellow |
| <u>_exptl_crystal_size_max</u> | 0.40 |
| <u>_exptl_crystal_size_mid</u> | 0.05 |
| <u>_exptl_crystal_size_min</u> | 0.05 |
| <u>_diffrn_ambient_temperature</u> | 100(2) |

5. ACKNOWLEDGEMENT

The crystallographer should be considered for co-authorship when the structural information is an important part of the paper and structural information has been derived mainly from the diffraction data. The crystallographer will help completing the X-ray part of the experiments and reviewing the manuscript. If structure determination was used only to confirm information obtained by other means (NMR, MS, etc.) and no structural details will be given in the paper, only acknowledgment is more appropriate, i.e. we acknowledge the Department of Chemistry X-ray Diffraction Facility and Dr. Chunhua Hu for his help with structure determination.